


COVID-19


Virus Inactivation



Viral RNA Preservation



Sampling & Transport Kit



NUCLISWAB®

Nasopharyngeal Swab + Collection Tube

NUCLISWAB® is a sampling and transport system used for collection and transport of viruses and cells from clinical samples or from environment, for isolation of nucleic acids that will be used in nucleic acid amplification tests.

NUCLISWAB® is one of the best adapted kits to detect **the coronavirus** meeting the recommendations for **COVID-19 Specimen Collection & Handling**: the virus inactivation due to the effect of the transport solution on the virus envelop and the preservation of the viral RNA due to the non effect on the viral capsid and a stable pH ensure ideal protection of nucleic acids.

 **SALUBRIS Inc.**

 www.salubrisinc.com

 Headquarters & U.S. Office
Phone (Toll Free): (800) 326-1484

 info@salubrisinc.com

 Europe Office:
Phone: +33 6 83 00 84 77